



Tires



SCP_137A



Tires

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1565361	Rim tape 22-559	1	--
2		1464199	PU-Tube (airless) 22"	1	--
3		1456300	Tube 22 x 1" Schrader valve	1	--
4		1456277	Tube 24 x 1" Schrader valve	1	--
5		1565708	Tube (airless) 24" x 1.3/8" for C63 tire	1	--
6		1592503	PU massiv ProSpeed, 24"	1	--
7		1593141	Tube with Schrader valve (23-559 - Marathon Plus, 28-559 - Right Run)	1	--
8		1592423	Tire Rightrun 25-489 black	1	--
9		1592428	Tire punction proof 22" x 1-3/8"	1	--
10		1565542	Tire SmartGuard Marathon Plus Evo 25-540	1	--
11		1591983	Tire punction proof 24" x 1-3/8"	1	--
12		1592199	Tire Rightrun 25-540 black	1	--
13		1593142	Tire Rightrun 25 x 1.1" (28-559)	1	--
14		1592429	Tire Marathon Plus 25-559 black	1	--