


Pushrim



Pushrim

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1520684	Screw e-motion M15	1	--
2		1592156	Cover for Sensor	1	--
3		1592161	Support for Spring	1	--
4		1592164	Sensor mounted	1	--
5		1592215	Pushrim 24" Stainless steel	1	--
6		1592216	Pushrim 24" Stainless steel with coating	1	--
7		1592217	Pushrim 24" Carbolife Curve L	1	--
8		1592218	Pushrim 24"" Carbolife Quadro	1	--
9		1592246	Pushrim 22" Stainless steel	1	--
10		1592247	Pushrim 22" Stainless steel coated	1	--
11		1592253	Pushrim 25" stainless steel	1	--
12		1592254	Pushrim 25" Carbolive Curve L	1	--
13		1592416	Spring	1	--
14		1592417	Cover for bracket	1	--
15		1592934	Sliding washer for holder of pushrim	1	--
16		1593234	Screw DIN7984 M4x5	1	--
17		1520680	Bushing	1	--