


# Anti-Tippers with Jack-Up Function since 2020



SCP\_141A

## Anti-Tippers with Jack-Up Function since 2020

| Item Number | Note | Part Number | Description   | Quantity Per |          |
|-------------|------|-------------|---|--------------|----------|
|             |      |             |   | Package      | Assembly |
| 1           |      | 1456183     | Washer DIN 125 - A 6,4  | 1            | --       |
| 2           |      | 1456233     | Cylinder head screw M 6 x 20  | 1            | --       |
| 3           |      | 1456406     | Quickpin quick release bolt   | 1            | --       |
| 4           |      | 1459636     | Hexagon nut ISO 4032 - M8 vz  | 1            | --       |
| 5           |      | 1459944     | Washer DIN 988 - 10 x 16 x 1  | 1            | --       |
| 6           |      | 1463698     | DIN 988 - 10 x 16 x 0,3 - St  | 1            | --       |
| 7           |      | 1467811     | Friction bearing Igus form F GFM-1212-035                                 | 1            | --       |
| 8           |      | 1467845     | Support angle right for anti-tippers                                      | 1            | --       |
| 9           |      | 1467846     | Support angle left for anti-tippers                                       | 1            | --       |
| 10          |      | 1467943     | Threaded pin M 8 x 45 DIN 913 - 45 H                                      | 1            | --       |
| 11          |      | 1489219     | Hardware set A for anti-tippers small version 2 x 15mm shoe incl. screws  | 1            | --       |
| 12          |      | 1489220     | Hardware set for anti-tippers expanded version 4 x 15mm shoe incl. screws | 1            | --       |
| 13          |      | 1489221     | Hardware set for anti-tippers expanded version 4 x 12mm shoe inc. screws  | 1            | --       |
| 14          |      | 1489365     | Wheel for anti-tippers complete   | 1            | --       |
| 15          |      | 1489366     | Slide rod anti-tippers mounted  | 1            | --       |
| 16          |      | 1519603     | Securing disk DIN 6799-Gr.8 vz  | 1            | --       |
| 17          |      | 1593206     | Anti-tippers with jack-up function (Pair) black                           | 2            | --       |