



Battery Charger and Cables



SHEP_241A

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60102598	Kit, Charger, Battery, 8 Amp. 115vac (Off Board) Single Mode Round Plug w/ AC Power Line Cord	1	1
2		60102582	Kit, Cord, AC Power, 115VAC, 3 Prong On-Board (6' 6-1/2" L)	1	1
3	B	60102824	Kit, Cable, Battery Extension	1	1

NOTE: A - CTE charger is standard on-chair. Used on Gel Cell or AGM (Absorbed Glass Mat) batteries.
 B - Used on GP24 Post Terminals and (22NF Post Terminals Before 02/14/22).