


Cables - Bus (PWH), Adapter, Sensor, Actuator



SHEP_347B

Cables - Bus (PWH), Adapter, Sensor, Actuator

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60124878	Kit, Cable, Accessory w/ Plug	1	1
			Cable, Accessory w/ Plug 420mm 16-17/32"	1	1
			Sheet, Insert Parts Service	1	1
2		60124378	Harness, Wire, TYCO-ACT 270mm 10-5/8"	1	1
			Harness, Wire, TYCO-ACT 270mm 10-5/8"	1	1
			Sheet, Insert Parts Service	1	1
2		60102526	Kit, Adapter Cable Assembly, TYCO-ACT 500mm 19-11/16"	1	1
			Adapter Cable Assembly, TYCO-ACT 500mm 19-11/16"	1	1
			Sheet, Insert Parts Service	1	1
3a		60102546	Kit, Adapter Cable Assembly, DX2-ACT4 To Molex 2 Pin	1	1
			Adapter Cable Assembly, DX2-ACT4 To Molex 2 Pin	1	1
			Sheet, Insert Parts Service	1	1
3a	C	60102545	Kit, Actuator Cable, LNX, (MPELR - Left)	1	1
			Actuator Cable, LNX, (MPELR - Left)	1	1
			Sheet, Insert Parts Service	1	1
3a	C	60105177	Kit, Actuator Cable (MPELR - Right)	1	1
			Actuator Cable (MPELR - Right)	1	1
			Sheet, Insert Parts Service	1	1
3b	D	60132729	Kit, Harness Wire Extension, 650mm 25-19/32"	1	1
			Harness Wire Extension, 650mm 25-19/32"	1	1
			Sheet, Insert Parts Service	1	1
4		60102528	Kit, LiNX Bus Cable (PWH) 300mm 11-13/16" (GLM-BUS030)	1	1
			LiNX Bus Cable (PWH) 300mm 11-13/16" (GLM-BUS030)	1	1
			Sheet, Insert Parts Service	1	1
4		60102529	Kit, LiNX Bus Cable (PWH) 700mm 27-9/16" (GLM-BUS070)	1	1
			LiNX Bus Cable (PWH) 700mm 27-9/16" (GLM-BUS070)	1	1
			Sheet, Insert Parts Service	1	1
4		60102536	Kit, LiNX Bus Cable (PWH) 1200mm 47-1/4" (GLM-BUS120)	1	1
			LiNX Bus Cable (PWH) 1200mm 47-1/4" (GLM-BUS120)	1	1
			Sheet, Insert Parts Service	1	1
4		60102530	Kit, LiNX Bus Cable (PWH) 1500mm 59-1/16" (GLM-BUS150)	1	1
			LiNX Bus Cable (PWH) 1500mm 59-1/16" (GLM-BUS150)	1	1
			Sheet, Insert Parts Service	1	1
4		60102531	Kit, LiNX Bus Cable (PWH) 1700mm 66-59/64" (GLM-BUS170)	1	1
			LiNX Bus Cable (PWH) 1700mm 66-59/64" (GLM-BUS170)	1	1
			Sheet, Insert Parts Service	1	1

Cables - Bus (PWH), Adapter, Sensor, Actuator

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
4		60124987	Kit, LiNX Bus Cable (PWH) 1000mm 39-3/8" (GLM-BUS100)	1	1
			LiNX Bus Cable (PWH) 1000mm 39-3/8" (GLM-BUS100)	1	1
			Sheet, Insert Parts Service	1	1
4		60126411	Kit, LiNX Bus Cable (PWH) 500mm 19-11/16" (GLM-BUS 050)	1	1
			LiNX Bus Cable (PWH) 500mm 19-11/16" (GLM-BUS 050)	1	1
			Sheet, Insert Parts Service	1	1
5a		60124379	Kit, Harness, Wire Extension 1040mm 40-15/16"	1	1
			Harness, Wire Extension 1040mm 40-15/16"	1	1
			Sheet, Insert Parts Service	1	1
5a	C	60124384	Harness, Wire Extension 650mm 25-9/32"	1	1
			Harness, Wire Extension 650mm 25-9/32"	1	1
			Sheet, Insert Parts Service	1	1
5b	D	60132730	Kit, Harness Wire Extension, 1473mm 57-63-64"	1	1
			Harness Wire Extension, 1473mm 57-63-64"	1	1
			Sheet, Insert Parts Service	1	1
6	A	60124849	Kit, C629 Junction Box Extension Cable, Long 1010mm 39-49/64"	1	1
			C629 Junction Box Extension Cable, Long 1010mm 39-49/64"	1	1
			Sheet, Insert Parts Service	1	1
7	B	60124879	Kit, C629 Junction Box Extension Cable 400mm 15-3/4"	1	1
			C629 Junction Box Extension Cable 400mm 15-3/4"	1	1
			Sheet, Insert Parts Service	1	1
8		60126412	Kit, LiNX Bus Extension Cable 640mm 25-13/64" (GLM-BUSEXT 064)	1	1
			LiNX Bus Extension Cable 640mm 25-13/64" (GLM-BUSEXT 064)	1	1
			Sheet, Insert Parts Service	1	1
NOTE: A - Long Extension Cable used with Elevate Option. B - Extension Cable used with Non-Elevate Option. C - Used before 01/24/22 D - Used after 01/23/22 NOTE: Refer to the Wiring Diagrams for detail application and configurations of Cables and Harnesses located at www.invacare.com/homecare under Product Literature Search/ Service Manual. 60108795 "LINX Control System, AVIVA FX", 60107431 "AVIVA FX Base", 60108797 "Ultra Low Maxx by Motion Concepts Modular Positioning System".					