


# C629M Junction Box, Shroud Cover and Hardware



## C629M Junction Box, Shroud Cover and Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60132727	Kit, Junction Box, C629M	1	1
2	B	60102487	Kit, Seat Frame Shroud Assembly	1	1
3	C	60132738	Kit, C629M Junction Box Extension Cable, Long 1010mm 39-49/64"	1	1
			C629M Junction Box Extension Cable, Long 1010mm 39-49/64"	1	1
			Sheet, Insert Parts Service	1	1
4	D	60132739	Kit, C629M Junction Box Extension Cable, 400mm 15-3/4"	1	1
			C629M Junction Box Extension Cable, 400mm 15-3/4"	1	1
			Sheet, Insert Parts Service	1	1

NOTE: A - C629 Junction Box includes the Hook and Loop Adhesive Strips.  
 B - Use existing hardware when attaching the Seat Frame Shroud.  
 C - Long Extension Cable used with Elevate Option.  
 D - Extension Cable used with Non-Elevate Option.