

Accessories Ocean (24") Ergo (XL) Before 02/01/21



SPC_063A

Accessories Ocean (24") Ergo (XL) Before 02/01/21

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		AP1603155	Splash guard - Ocean Ergo Series	1	1
2	A	AP1603142	Standard soft seat - Ocean Ergo Series	1	1
3		AP1603143	Soft seat insert - Ocean Ergo Series	1	1
4	B	AP1603145	Small soft seat - Ocean Ergo Series	1	1
5		AP1603152	Small soft seat insert - Ocean Ergo Series	1	1
6		1470081	Chest Belt Padded w/ Stainless Hardware	1	1
6		1470082	Pelvis Belt Padded w/Stainless Hardware	1	1
7	C	1470077	Armrest lock (pair)	2	1
8		1535077	Back Lateral (EA)	1	1
9		AP1643543	Legrest Cpl., right - Ocean / Ocean Ergo Series	1	1
9		AP1643572	Legrest Cpl., left - Ocean / Ocean Ergo Series	1	1
10		AP1643573	Thigh support, right - Ocean / Ocean Ergo Series	1	1
10		AP1643574	Thigh support, left - Ocean / Ocean Ergo Series	1	1
11		AP1604136	Comfort soft backrest - Ocean (24") Ergo	1	1
11		AP1604016	Comfort soft backrest - Ocean (24") Ergo XL	1	1
12		1558769	Soft cushion w/cover blue - Ocean / Ocean Ergo Series	1	1
13		AP1602078	Collection Pan holder - Ocean Ergo Series	1	1
14		10230	Collection Pan with Lid	1	1
15		AP1601986	Front anti-tipper - Ocean (24") Ergo (XL)	1	1
16		16379	Self Propelling Wheels Kit complete including brakes	1	1
17		AP1604015	XL backrest kit - Ocean (24") Ergo	1	1
18		AP1606985	Backrest tube only -Ocean VIP Ergo (separate order needed for armrest, backrest cover, headrest)	1	1

NOTE: A - Center key shaped oval cut out is used for standard adult size individual. 11-13/16" is the approximate widest width of the oval cut out.
 B - Center key shaped oval cut out is used for smaller and young adults. 9-13/16" is the approximate widest width of the oval cut out.
 C - Only for Ocean Ergo with XL-armrests (1470086)