

Accessories OceanErgoVIP After 01/31/21



Accessories OceanErgoVIP After 01/31/21

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		AP1603155	Splash guard - Ocean Ergo Series	1	1
2	A	AP1603142	Standard soft seat - Ocean Ergo Series	1	1
3		AP1603143	Soft seat insert - Ocean Ergo Series	1	1
4	B	AP1603145	Small soft seat - Ocean Ergo Series	1	1
5		AP1603152	Small soft seat insert - Ocean Ergo Series	1	1
6		AP1603154	Universal soft seat - Ocean Ergo Series	1	1
7		AP60121311	Special Soft Seat, Ocean Ergo Series	1	1
8		AP1605002	Protection Cvr, Variable soft seat	1	1
9	D	1544469	Variable soft seat cpl. - Ocean / VIP / Ergo Series	1	1
9		1545659	Clip, Variable Softseat	1	2
10		AP60123732	Safety Bar 2.0 (for armrest 2.0)	1	1
10		1470078	Safety Support Bar , curved	1	1
11		AP1643760	Safety Lock, 2.0 (1 lock) (for armrest 2.0)	1	2
12		1470081	Chest Belt Padded w/ Stainless Hardware	1	1
12		1470082	Pelvis Belt Padded w/Stainless Hardware	1	1
13		1535077	Back Lateral (EA)	1	2
14	E	1470077	Armrest lock (pair)	2	1
15		1470086	Armrests, Extra Large (Pair)	2	1
16		AP1643543	Legrest Cpl., right - Ocean / Ocean Ergo Series	1	1
16		AP1643572	Legrest Cpl., left - Ocean / Ocean Ergo Series	1	1
17		AP1643573	Thigh support, right - Ocean / Ocean Ergo Series	1	1
17		AP1643574	Thigh support, left - Ocean / Ocean Ergo Series	1	1
18		AP1603606	Comfort soft backrest - Ocean (Dual-)VIP Ergo	1	1
19		1575299	Ocean Contour Hygienic Fixed Back Support	1	1
20		1558769	Soft cushion w/cover blue - Ocean / Ocean Ergo Series	1	1
21		AP1606976	Reclining backrest - OceanErgoVIP	1	1
22		1531465	Ergonomic headrest pad, small	1	1
23		1531608	Ergonomic headrest pad, large	1	1
24		1531613	Adjustable headrest Bracket	1	1
25		AP1643748	Armrest Pad, 2.0 (1 Pad)	1	2
26		AP1602078	Collection Pan holder - Ocean Ergo Series	1	1
27		10230	Collection Pan with Lid	1	1
28		AP1603992	Front and back anti-tipper - Ocean (Dual-)Vip Ergo	1	1
29	C	AP60129156	Tilt Stopper Screw Set	1	1

NOTE: A - Center key shaped oval cut out is used for standard adult size individual. 11-13/16" is the approximate widest width of the oval cut out.
 B - Center key shaped oval cut out is used for smaller and young adults. 9-13/16" is the approximate widest width of the oval cut out.
 C - Tilt Stopper implementation Serial #21DVP300001 for models w/o Caster Directional Lock and 21DVP310001 w/ Caster Directional Lock.
 D - It is necessary to use item 8 with Item 9. Item 8 is used to avoid clamping risk to the fingers when using item 9.
 E - Only for Ocean Ergo with XL-armrests (1470086)

Item numbers not listed are non-saleable.