


# Spare parts Ocean Ergo (XL) After 01/31/21

Item numbers not listed are non-saleable.



# Spare parts Ocean Ergo (XL) After 01/31/21

Item numbers not listed are non-saleable.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		SP1607166	Backrest w/o arm. / Hdwr. - Ocean (24") Ergo	1	1
2	A,B	SP1647584	Armrest Assembly 2.0 w/ Hardware, Left	1	1
3		SP1647608	Screws (Set) for Armrest 2.0	1	2
4		SP1607167	Backrest XL w/o arm. / Hdwr. - Ocean (24") Ergo XL Armrest	1	1
5		SP1643752	Cover, Armrest 2.0	1	2
6		SP1607743	Seat plate - Ocean Ergo Series	1	1
7		SP1607371	Side Frame w / 5" Casters - Ocean (24") Ergo (XL)	1	1
8		16334	Heel Loop (Pair)	2	1
9		SP1607170	Hdwr., heel strap - Ocean (24") Ergo (XL)	1	2
10		SP1607168	Base stabilizer bar w / Hdwr. - Ocean (24") Ergo XL Footrest cpl.,	1	1
11		SP1643541	Footrest cpl., Left - Ocean (XL) / (Dual-)VIP / Ocean Ergo Series	1	1
13		1470689	Footrest Plate	1	2
14		SP1643542	Footrest cpl., Right - Ocean (XL) / VIP / Ocean Ergo Series	1	1
15		SP1607628	Caster 5" - Ocean Ergo Series	1	4
16		SP1607546	Side Frame w / 5" Dir. Caster - Ocean (24") Ergo (XL)	1	1
17		SP1607742	Clip, footrest - Ocean Ergo Series	1	2
18		SP1607741	Caster 5", dir. lock - Ocean Ergo Series	1	1
19		SP1607169	Plugs, Side Frame - Ocean (24") Ergo (XL)	1	1
20		16256	Side Frame Adj. Clip	1	4
21		1542900	Bkrst Hdwr. - Ergo / Ocean (Dual-)VIP / XL	1	1
22	A,B	SP1647583	Armrest Assembly 2.0 w/ Hardware, Right	1	1
23		AP1603031	Backrest Cover - Ocean Ergo/ SP	1	1
24		AP1603577	Backrest Cover - Ocean Ergo XL/ SPXL	1	1
25	C,D	SP60129114	Seat Height Lowering Screw Set	1	1

NOTE: A - Armrest 2.0 implementation Serial #21BVP100229 for Ergo model w/o Caster Directional Lock and 21EVP110000 w/ Caster Directional Lock.  
B - Armrest 2.0 implementation Serial #21BVP120021 for ErgoXL model w/o Caster Directional Lock and 21EVP130001 w/ Caster Directional Lock.  
C - Seat Height Lowering Screw Kit implementation Serial #21FVP100001 for Ergo model w/o Caster Directional Lock and 21EVP110000 w/ Caster Directional Lock.  
D - Seat Height Lowering Screw Kit implementation Serial #21FVP120001 for ErgoXL model w/o Caster Directional Lock and 21FVP130001 w/ Caster Directional Lock.