


Seat Extended Legrest

88



Seat Extended Legrest 88

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	04529X035-NLA	NLA - Pad, Legrest (12" W)	1	1
1	A	04529X036-NLA	NLA - Pad, Legrest (14" W)	1	1
1	A	04529X037-NLA	NLA - Pad, Legrest (16" W)	1	1
2		1027095	Plug Button, Black (7/8")	1	1
3		48127X001-NLA	NLA - Screw, Phillips Round Head (#10-32 x 1/2")	1	1
4		01058X005-NLA	NLA - Tube, Seat Extension	1	1
5		40188X009	Button, Single Snap	1	1
6	A	04706X082-NLA	NLA - Seat Extended Legrest Assembly (12" W)	1	1
6	A	04706X072-NLA	NLA - Seat Extended Legrest Assembly (14" W)	1	1
6	A	04706X074-NLA	NLA - Seat Extended Legrest Assembly (16" W)	1	1
NOTE: A - Upholstery color is only available in U16 Light Blue Heavyweight Vinyl.					