


Solid Folding Seat



Solid Folding Seat

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1016766	Bracket, Hinge, Chrome	1	2
1		1020225-NLA	NLA - Bracket, Hinge, Chrome (14" W Chair)	1	2
2		1013218	Washer (#10 Flat)	1	8
3		7930360-NLA	NLA - Package, Screw, Phillips Pan Head (#10-32 x 11/16")	25	8
4	A,B	8881020216-NL A	NLA - Upholstery, Solid Seat	1	1
4	A,C	8881020220-NL A	NLA - Upholstery, Solid Seat	1	1
4	A,D	8881020227-NL A	NLA - Upholstery, Solid Seat	1	1
5		1016767-NLA	NLA - Hook Bracket	1	2
5		1020224	Hook Bracket, Hook (14" W)	1	2
6	A,B	8881036890-NL A	NLA - Solid Folding Seat Assembly	1	1
6	A,C	8881036891-NL A	NLA - Solid Folding Seat Assembly	1	1
6	A,D	8881036892-NL A	NLA - Solid Folding Seat Assembly	1	1
NOTE: A - Specify upholstery color, see page titled "Packaging and Upholstery Colors" B - 12" Wide x 11" Deep C - 14" Wide x 12" Deep D - 16" Wide x 14" Deep					