


Back Canes and Hardware



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1006229	Grip, Foam (8")	1	1
2		1042351	Bushing (3/4")	1	1
3		1042401	Tube, Straight	1	1
3		1020545-NLA	NLA - Cane, Back with Bend	1	1
4		1008432	Screw, Socket Button Head (1/4-20 x 1-1/4")	1	2
5		7930364-NLA	NLA - Package, Locknut (1/4-20)	25	2
6	A	1028372	Back Cane Assembly, 10 Degree Bend	1	1
6	A	1042402-NLA	NLA - Back Cane Assembly, Straight	1	1

NOTE: A - Includes: Items 1-3