



Rear Wheels (without Handrims)



7 Spoke Composite Wheel



9 Spoke Composite Wheel

Rear Wheels (without Handrims)

Wheel Style	Wheel Size	Tire Type	Single Wheel Assembly	Tire	Tube Insert	Bearing	Bearing Spacer
U1763BK See Note(s) C	20x1-1/4"	Urethane	1038555 See Note(s) C, B, A	1037774	N/A	40117X024	1015756
1763BK	20x1-3/8"	Pneumatic	1034126 See Note(s) B, A	1017519	1017520	40117X024	1015756
4763BK	20x1-3/8"	Pneumatic w/ Flat Free Insert	1034127 See Note(s) C, B, A	1017519	1023322	40117X024	1015756
U1764BK	22x1-1/4"	Urethane	1032115 See Note(s) C, B, A	1028250	N/A	40117X024	1015756
1764BK	22x1-3/8"	Pneumatic	1027149 See Note(s) B, A	1018544	1018545	40117X024	1015756
4764BK	22x1-3/8"	Pneumatic w/ Flat Free Insert	1027150 See Note(s) A, B	1018544	1023323	40117X024	1015756
U2244BK 7-Spoke	24x1-1/4"	Urethane	1032117 See Note(s) A, C	1026904	N/A	40117X024	1015756
2244BK	24x1-3/8"	Pneumatic	1025123 See Note(s) A	40176X060	40176X061	40117X024	1015756
4244BK	24x1-3/8"	Pneumatic w/ Flat Free Insert	1025124 See Note(s) A, C	40176X060	1010906	40117X024	1015756
U2222 9-Spoke	24 x 1"	Urethane	1044193 See Note(s) C, A, D	1039304	N/A	40117X024	1015756
U2222 9-Spoke	24 x 1"	Urethane	1080378 See Note(s) F, E, A, C	1080400 See Note(s) F	N/A	40117X024	1015756
U2222C 9-Spoke	24 x 1"	Urethane Corded	1122726 See Note(s) A	1118209 See Note(s) G	N/A	40117X024	1015756

NOTE: A - Includes tire, tube/insert, bearing, and spacer; Does not include handrim.
 B - When ordering 20" or 22" wheels, a wheel lock shoe spacer may also be required. See page titled "Wheel Locks".
 C - Composite Wheels with Flat Free Inserts or Urethane Tires have a plug which caps an opening in the rim. See page titled "Rear Wheel Hardware" for plug, part number 1030936.
 D - Before 11/20/98
 E - After 11/19/98
 F - Non-Corded Tire replaced with a Corded Tire as of 11/16/09. Refer to Wheel Style U2222C part #1122726 for a wheel assembly or part #1118209 for a tire. For chairs manufactured before 11/16/09, it is recommended to order (2) each in quantity to match.
 G - Corded Tires maybe more difficult to stretch over the rim. The complete wheel assembly that has the tire and wheel rim preassembled maybe advisable.