




Tash and Therafin Accessories




P5820 - Red Soft Switch
P5821 - Green Soft Switch




✦ P5018 - Wafer Board




✦ P5026 - Star Control



✦ P5020 - Mini Joystick with Push



P4337 - Multi-Ability Cable



P4353 - CA5 Multi-Connector
Box



P5970 - "U" Clamp

✦ Must be used with display w/ switches (1108489)

Tash and Therafin Accessories

Therafin

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
	B	1142362	Bib Assembly (PKG32669)	1	1
	B	P21456-NLA	NLA - Adaptor Clamp - Tash	1	1
	B	P21457-NLA	NLA - Adaptor Clamp - Invacare	1	1
	B	P21460-NLA	NLA - Midline Swingway System	1	1
	B	1140824	Kit, Sip-N-Puff, 16" Fixed (PKG32666) Includes: (3) Air Tube Straws, (2) Saliva Traps, Tubing, Pivot Shaft/Flexible Metal Gooseneck Assembly, and Extruded Figure "8" Clamp to fit 7/8" Wheelchair Frame Tubing. Extruded Back Cane Adapter Clamp.	1	1
	A	P5820-NLA	NLA - Switch, Red Soft	1	1
	A	P5821	NLA - Switch, Green Soft	1	1
	A	P5018-NLA	NLA - Wafer Board (5018)	1	1
	A	P5026-NLA	NLA - Star Control	1	1
	A	P5020-NLA	NLA - Mini Joystick with Push (5020)	1	1
	A	P4337-NLA	NLA - Multi-Ability Cable	1	1
	A	P4353-NLA	NLA - CA5 Multi-Connector Box	1	1
	A	P5970-NLA	NLA - Clamp, "U"	1	1
	B	1108489	Display with Switches	1	1
NOTE: A - See Pictorial B - Not Shown					
NOTE: For more information on Tash and Therafin products, please contact one of our Rehab Specialists.					

Tash and Therafin Accessories

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		P5820-NLA	NLA - Switch, Red Soft	1	1
1		P5821	NLA - Switch, Green Soft	1	1
2	C	P5018-NLA	NLA - Wafer Board (5018)	1	1
3	C	P5026-NLA	NLA - Star Control	1	1
4	C	P5020-NLA	NLA - Mini Joystick with Push (5020)	1	1
5		P4337-NLA	NLA - Multi-Ability Cable	1	1
6		P4353-NLA	NLA - CA5 Multi-Connector Box	1	1
7		P5970-NLA	NLA - Clamp, "U"	1	1
	B	1142362	Bib Assembly (PKG32669)	1	1
	B	P21456-NLA	NLA - Adaptor Clamp - Tash	1	1
	B	P21457-NLA	NLA - Adaptor Clamp - Invacare	1	1
	B	P21460-NLA	NLA - Midline Swingway System	1	1
	A,B	1140824	Kit, Sip-N-Puff, 16" Fixed (PKG32666) Includes: (3) Air Tube Straws, (2) Saliva Traps, Tubing, Pivot Shaft/Flexible Metal Gooseneck Assembly, and Extruded Figure "8" Clamp to fit 7/8" Wheelchair Frame Tubing. Extruded Back Cane Adapter Clamp.	1	1
	B	1108489	Display with Switches	1	1

NOTE: A - See items included.
B - Not Shown
C - Must be used with display with switches part #1108489.