



# Wheel Locks with Hill Holder



Attaching hardware

## Wheel Locks with Hill Holder

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		7930364-NLA	NLA - Package, Locknut (1/4-20)	25	1
2		7930367-NLA	NLA - Locknut (1/4-28)	25	3
3		7930357-NLA	NLA - Washer, Nylon (1/4 x 1/2 x 1/32")	25	7
4		1015445-NLA	NLA - Bracket, Mounting, Chrome - Right	1	1
4		1015446-NLA	NLA - Bracket, Mounting, Chrome - Left	1	1
5		40157X004	Spacer, Nylon (1/4 x 7/16 x 1/8")	1	1
6		60274X012-NLA	NLA - Screw, Hex Head (1/4-28 x 3/4")	1	3
7		1017861-NLA	NLA - Link, Chrome	1	1
8		47376X024	Tip, Rubber Wheel Lock Handle, Black	1	1
9	E	1016588-NLA	NLA - Handle, Push, Adult, Chrome - Right	1	1
9	E	1016589-NLA	NLA - Handle, Push, Adult, Chrome - Left	1	1
9	E	1016590-NLA	NLA - Handle, Pull, Hemi, Chrome - Right	1	1
9	E	1016591-NLA	NLA - Handle, Pull, Hemi, Chrome - Left	1	1
10		1010398	Washer, Nylon (1/4 x 1/2 x 1/16")	1	3
11		1007017-NLA	NLA - Screw, Hex Head (1/4-20 x 1-1/4")	1	1
12	H	1016892-NLA	NLA - Hill Holder Assembly - Right	1	1
12	H	1016893-NLA	NLA - Hill Holder Assembly - Left	1	1
13	A,B,C	1034212	Wheel Lock Assembly, Hill Holder, Push to Lock - Right	1	1
13	A,B,C	1034213	Wheel Lock Assembly, Hill Holder, Push to Lock - Left	1	1
13	A,B,D	1034214	Wheel Lock Assembly, Hill Holder, Pull to Lock - Right	1	1
13	A,B,D	1034215	Wheel Lock Assembly, Hill Holder, Pull to Lock - Left	1	1
14	E	7930351-NLA	NLA - Nut, Hex Flange (1/4-20)	25	1
14	F	1021403	Nut, Hex Flang (1/4-20)	1	1
15	E	7930309-NLA	NLA - Bolt, Carriage (1/4-20 x 1-1/2" )	25	1
15	F	1021402	Bolt, Carriage (1/4-20 x 1-1/4")	1	1

NOTE: A - Hill holder Wheel Locks not for use with Urethane tires.  
 B - Includes items 1 thru 12  
 C - Push to Lock - Adult  
 D - Pull to Lock - Hemi  
 E - Chrome Chairs  
 F - Black or Color Chairs  
 H - Part #1016892 and 1016893 are no longer available as a service part individually. Order item #13 as a replacement.