



Deluxe Telescoping Bed Rail and Crossbrace Assembly 6625/6625L

6625 was discontinued 2/4/03.
6625L was discontinued 9/21/94.

Model 6625



Model 6625L (Discontinued) (Fits Wider Twin Beds-Deck Width 33-1/4" - 36-1/2")



Deluxe Telescoping Bed Rail and Crossbrace Assembly 6625/6625L

**6625 was discontinued 2/4/03.
6625L was discontinued 9/21/94.**

| Item Number | Note | Part Number | Description | Quantity Per | |
|---|------|---------------|---|--------------|----------|
| | | | | Package | Assembly |
| 1 | | 16320X006-NLA | NLA - Rail, Outside | 1 | 2 |
| 1 | | 16320X007-NLA | NLA - Rail, Inside | 1 | 2 |
| 2 | | 1057287 | Plug Button, Grey (1") | 1 | 8 |
| 3 | | 40314X001-NLA | NLA - Ring, Retaining | 1 | 4 |
| 4 | A,D | 16330X000-NLA | NLA - Crossbrace Assembly, Grey | 1 | 2 |
| 4 | E | 1010422-NLA | NLA - Crossbrace Assembly | 1 | 2 |
| 5 | D | 1007727 | Channel Only, Granite | 1 | 2 |
| 5 | E | 1007729-NLA | NLA - Channel Only | 1 | 2 |
| 9 | | 60290X002 | Locknut (5/16-18) | 1 | 4 |
| 10 | D | 16308X000 | Washer, Flat (5/16 x 11/16 x 1/16") | 1 | 4 |
| 10 | E | 16347X000-NLA | NLA - Washer, Nylon (3/16 x 7/16 x 1/8") | 1 | 4 |
| 11 | B | 16326X000 | Barrel Assembly, Left Head End / Right Foot End (Includes: Plunger Knob and Pin) | 1 | 2 |
| 11 | B | 16356X000 | Barrel Assembly, Right Head End / Left Foot End (Includes: Plunger Knob and Pin) | 1 | 2 |
| 12 | | 19000X023 | Spring, Plunger (23/64 x 5/8") | 1 | 4 |
| 13 | | 19000X012 | Pin, Plunger | 1 | 4 |
| 14 | | 40014X020-NLA | NLA - Knob, Plunger | 1 | 4 |
| 15 | C | 1028789-NLA | NLA - Kit, Frame Bar (Pair) | 1 | 1 |
| 6 | | | Bolt, Carriage (1/4-20 x 1-1/4") | 1 | 2 |
| 7 | | | NLA - Clamp Bar | 1 | 2 |
| 8 | | | Nut, Wing (1/4-20) | 1 | 2 |
| NOTE: A - Includes Items 5-11 B - Includes Items 3, 12-14 C - Includes Items 6-8 D - Used on 6625 E - Used on 6625L | | | | | |