


# Oxygen Holder



# Oxygen Holder

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1108517-NLA	NLA - Bracket, Oxygen Bottle Mounting	1	1
2		1044939-NLA	NLA - Screw, Hex Head (5/16-18 x 1-3/4")	1	4
3		1017470	Clamp, Tube (1")	1	2
4		1028769	Locknut (5/16-18)	1	4
5		1108550-NLA	NLA - Oxygen Holder Assembly w/ Instruction Sheet	1	1
	A,B	1110522-NLA	NLA - Instruction Sheet	1	1

NOTE: A - Not Shown  
 B - Item is no longer available from the Invacare Service Parts Department. Refer to the Invacare Website under "Product Literature Search". Enter the part number in the "Keywords" field and download a free PDF version copy.