


Seat Pans



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1008495	Screw, Socket Button Head (1/4-20 x 1-1/2")	1	4
		1025195	Locknut (1/4-20)	1	4
2		See Chart	#2 Seat Pans	1	1

NOTE: A - Not Shown

#2 Seat Pans

Seat Depth	10" Wide	11" Wide	12" Wide	13" Wide	14" Wide	15" Wide	16" Wide
10-13"	1071198	1071199	1073589	1073590	1073591	1073592	1073593

Seat Pans

#2 Seat Pans

Seat Depth	10" Wide	11" Wide	12" Wide	13" Wide	14" Wide	15" Wide	16" Wide
14-16"	1073594	1073595	1073596-NLA	1073597	1073598	1073599-NLA	1073600