





Side Frames and Hardware (Before 10/1/00)



Side Frames and Hardware (Before 10/1/00)

Item #1 - Side Frames and Hardware

Frame Type	Non-Recliner	Recliner
Adult - Right	8881058593-NLA See Note(s) A, B	8881062457-NLA See Note(s) A, B
Adult - Left	8881058594-NLA See Note(s) A, B	8881062458-NLA See Note(s) A, B
Hemi - Right	8881058599-NLA See Note(s) B, A	8881062459-NLA See Note(s) B, A
Hemi - Left	8881058600-NLA See Note(s) A, B	8881062460-NLA See Note(s) A, B

NOTE: A - Specify frame finish when ordering. See page titled "Frame Finishes, Touch-Up Paints and Upholstery Colors."
B - Decals not included with frame. See page titled "Packaging, Literature and Decals."

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	Item #1 - Side Frames and Hardware	1	1
2		1090768	Kit, Lock Handle w/Rivets (Pair)	2	1
3		47471X074	Spring, Arm Lock	1	1
4		1028594	Pin, Arm Lock, Black	1	1
5		40117X006	Bearing (5/8")	1	2
6		1059371-NLA	NLA - Cap, End, Oval	1	1
7	A	1085124-NLA	NLA - Screw, Socket Flat Head (#10-32 x 1-3/4")	1	1
8		1007541	Clamp, Seat	1	2
9		1003545	Screw, Socket Head (1/4-20 x 1-1/4")	1	2
10	A	1085134-NLA	NLA - Screw, Socket Flat Head (#10-32 x 1-5/8")	1	1
11		1060099-NLA	NLA - Hanger, Front Rigging, Black - Right	1	1
11		1060100-NLA	NLA - Hanger, Front Rigging, Black - Left	1	1
12		1027097	Plug Button, Black (3/4")	1	1
13		1025736	Locknut (1/4-20)	1	2
14		1025731	Locknut (#10-32)	1	2
15		1058520-NLA	NLA - Fastener, Panel Snap-In	1	1

NOTE: A - Screws changed 9/16/99 and serve as a direct replacement for chairs manufactured before 9/16/99.