


# Back Canes and Attaching Hardware (After 9/30/00)



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1028372	Back Cane Assembly, 10 Degree Bend	1	2
1		1042402-NLA	NLA - Back Cane Assembly, Straight	1	2
2		1006229	Grip, Foam (8")	1	2
3		1071565-NLA	NLA - Tube, Rear Upright, Black	1	2
4		1026907	Washer (1/4 x 1/2 x 1/16")	1	4
5		1057955	Screw, Hex Head (1/4-20 x 1-1/4")	1	4
6		1071488-NLA	NLA - Insert, Threaded	1	2
7		1055348	Screw, Hex Head (1/4-28 x 1-3/4")	1	4
8		1026900	Washer (1/4 x 23/32 x 1/16")	1	8
9		1065362	Locknut (1/4-28)	1	1
10		1025195	Locknut (1/4-20)	1	1
11	A	1071564-NLA	NLA - Hardware, Back Cane Attaching (Pair)	2	1

NOTE: A - Includes Items 3-10