


Fixed Height, Removable Arms

Non-Reclining (Before 10/1/00 and All Recliners)



Fixed Height, Removable Arms Non-Reclining (Before 10/1/00 and All Recliners)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	8881015932-NL A	NLA - Pad, Arm, Vinyl, Full Length - Right (14" L x 2" W)	1	1
1	A	8881015933-NL A	NLA - Pad, Arm, Vinyl, Full Length - Left (14" L x 2" W)	1	1
1	A	8881011811-NL A	NLA - Pad, Arm, Vinyl, Desk Length (9" L x 2" W)	1	1
2		1055949	Package, Screw, Phillips Pan Head (#10-32 x 1-3/8")	2	1
3	B	8881057723	Desk Length Arm Tube (Color), Right	1	1
3	B	8881057724-NL A	NLA - Tube, Arm, Desk Length, Color - Left	1	1
3	C	1057725-NLA	NLA - Tube, Arm, Desk Length, Chrome - Right	1	1
3	C	1057726-NLA	NLA - Tube, Arm, Desk Length, Chrome - Left	1	1
3		1041891-NLA	NLA - Tube, Arm, Desk Length, Black - Right	1	1
3		1041895-NLA	NLA - Tube, Arm, Desk Length, Black - Left	1	1
3	B	8881057715-NL A	NLA - Tube, Arm, Full Length, Color - Right	1	1
3	B	8881057716-NL A	NLA - Tube, Arm, Full Length, Color - Left	1	1
3	C	1057717-NLA	NLA - Tube, Arm, Full Length, Chrome - Right	1	1
3	C	1057718-NLA	NLA - Tube, Arm, Full Length, Chrome - Left	1	1
3		1044513-NLA	NLA - Tube, Arm, Full Length, Standard, Black - Right	1	1
3		1044514-NLA	NLA - Tube, Arm, Full Length, Standard, Black - Left	1	1
4		1035472-NLA	NLA - Guard, Clothing, Desk Length	1	1
4		1044572-NLA	NLA - Guard, Clothing, Full Length (16 - 19" D)	1	1
5		1025729	Screw, Phillips Pan Head Tap (#8-32 x 3/8")	1	8
6	B	8881057729-NL A	NLA - Arm Assembly, Desk Length, Color - Right	1	1
6	B	8881057730-NL A	NLA - Arm Assembly, Desk Length, Color - Left	1	1
6	C	1057731-NLA	NLA - Arm, Desk Length without Armpad, Chrome - Right	1	1
6	C	1057732-NLA	NLA - Arm, Desk Length without Armpad, Chrome - Left	1	1
6		1057727-NLA	NLA - Arm Assembly, Desk Length, Chrome - Right	1	1
6		1057728-NLA	NLA - Arm Assembly, Desk Length, Chrome - Left	1	1
6	B	8881057719-NL A	NLA - Arm, Full Length without Armpad, Color - Right	1	1
6	B	8881057720-NL A	NLA - Arm, Full Length without Armpad, Color - Left	1	1
6	C	1057721-NLA	NLA - Arm, Full Length without Armpad, Chrome - Right	1	1
6	C	1057722-NLA	NLA - Arm, Full Length without Armpad (Chrome), Left	1	1
6		1044573-NLA	NLA - Arm Assembly, Full Length, Black - Right (Standard) (16-19" D)	1	1
6		1044574-NLA	NLA - Arm Assembly, Full Length, Black - Left (Standard) (16-19" D)	1	1

Fixed Height, Removable Arms Non-Reclining (Before 10/1/00 and All Recliners)

NOTE: A - Specify upholstery color when ordering. See page titled "Frame Finishes, Touch-Up Paints and Upholstery Colors."
B - Specify frame finish when ordering. See page titled "Frame Finishes, Touch-Up Paints and Upholstery Colors."
C - Chrome is NLA as of 3/22/99. Black replaced chrome. Order arm components in pairs to match.