


# P904A Legrest Supports



See ratchet and release handle assembly page.

**NOTE:** Calf pad and hardware not included with support asbly. See calf pad and hardware page.

## P904A Legrest Supports

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1027095	Plug Button, Black (7/8")	1	1
2		11919X000-NLA	NLA - Screw, Slotted Truss Head (3/8-16 x 1-1/2")	1	1
3		1058636-NLA	NLA - Washer, Nylon (25/64 x 7/8 x 1/32")	1	2
4		40363X000	Hex Nut (3/8"-16)	1	1
5		1025730	Locknut (5/16-24)	1	1
6		1026030	Bolt, Lug (5/16-24 x 1-3/8")	1	1
7		04526X040-NLA	NLA - Tube, Elevating Support, Chrome - Right	1	1
7		04526X041-NLA	NLA - Tube, Elevating Support, Chrome - Left	1	1
7	A	8881048851	Tube, Elevating Adjustable Support, Color - Right	1	1
7	A	8881048852	Tube, Elevating Adjustable Support, Color - Left	1	1
8		1027097	Plug Button, Black (3/4")	1	1
9		1058484-NLA	NLA - Tube, Legrest Support, Chrome - Right	1	1
9		1058485-NLA	NLA - Tube, Legrest Support, Chrome - Left	1	1
9	A	8881058664-NL A	NLA - Tube, Legrest Support, Color - Right	1	1
9	A	8881058665-NL A	NLA - Tube, Legrest Support, Color - Left	1	1
10		1022874	Bumper	1	1
11		1030979	Screw, Phillips Pan Head Tap (#10-32 x 3/4")	1	1
	A,B	8881058646-NL A	NLA - Legrest Support Assembly, Adult, Color - Right	1	1
	A,B	8881058647-NL A	NLA - Legrest Support Assembly, Adult, Color - Left	1	1
	B,C	1058488-NLA	NLA - Legrest Support Assembly, Adult, Silver Vein - Right	1	1
	B,C	1058489-NLA	NLA - Legrest Support Assembly, Adult, Silver Vein - Left	1	1

NOTE: A - Specify frame finish when ordering. See page titled "Frame Finishes, Touch-Up Paints and Upholstery Colors."  
 B - Not Shown  
 C - As of 7/6/07, chrome finish option is discontinued and no longer available.