

# Overbed Tables

If Models 6417 and 6418 have "Manufactured in China" on the column, please contact your Invacare Customer Service Representative at 1-800-333-6900 for Warranty Replacement information. For Models 6417 and 6418 "Not Manufactured in China" Please call 1-800-432-6322 and ask for Jeff Beatty (extension 242) or Kristie Corwin (extension 240).

Model 6416-W Auto Touch with Finger Tip Control (Discontinued)


Model 6417 Auto Touch (Shown)

Model 6418 Tilt Top

Model 6419 Heavy Duty



Model 6420 Vanity and Mirror  
(Discontinued)



# Overbed Tables

**If Models 6417 and 6418 have "Manufactured in China" on the column, please contact your Invacare Customer Service Representative at 1-800-333-6900 for Warranty Replacement information. For Models 6417 and 6418 "Not Manufactured in China" Please call 1-800-432-6322 and ask for Jeff Beatty (extension 242) or Kristie Corwin (extension 240).**

| Item Number   | Note  | Part Number | Description                  | Quantity Per |          |
|---|-------|-------------|------------------------------|--------------|----------|
|   |       |             |                              | Package      | Assembly |
|   | B,D   | 1134778-NLA | NLA - Package, Casters - OBT | 4            | 1        |
|   | A,C,D | 1144544     | Package, Casters             | 4            | 1        |
| NOTE: A - Includes (2) locking casters and (2) non-locking casters<br>B - Used with the 6417 and 6418, Before 10/06.<br>C - Used with 6417 and 6418, After 9/06.<br>D - The date code label is located on the column for the Overbed Table, when selecting castors for the 6417 and 6418 please note the date code in determining which caster is needed. |       |             |                              |              |          |