


# Seat Extension Kit

## 14" - 16" Deep Chairs Only



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1100972	Kit, Seat Extension (2")	1	1
2		1031095	Screws with Washers (#10-16 x 7/8")	4	2
3		1089278-NLA	NLA - Kit, Extension, Seat (2")	1	2
4		1100665	Tube, Seat Tension	1	2
5		1056328	Plug Button, Black (1")	1	2

**NOTE:** A - Seat extension kit includes items 2-4 but does not include a plug button.