


Rear Frame w/ Mounting Hardware



Rear Frame w/ Mounting Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1087632-NLA	NLA - Frame Assembly w/o Decal, Black - Rear	1	1
2		1071237	Plug Button, Square, Black (1")	1	4
3		1117326-NLA	NLA - Screw, Socket Head w/ Patch (5/16-18 x 7/8")	1	6
4		1069724	T-Nut (5/16-18)	1	7
5		1091869-NLA	NLA - Screw, Socket Flat Head w/ Patch (1/4-20 x 3/4")	1	4
6		1087652-NLA	NLA - Bar, Clamp, Color	1	2
7		1087732-NLA	NLA - Guide, Linear	1	2
8	A	1106550-NLA	NLA - Hardware, Rear Frame Mounting	1	1
9		1087602-NLA	NLA - Decal, Invacare Power Tilt Only	1	1
NOTE: A - Includes items 2-6; less Linear Guidance					