


Rear Wheels



7 Spoke Composite Wheel

NOTE: 22" Rear Wheels not available on 10" or 12" Deep non-reclining chairs.

Rear Wheels

7 Spoke Composite Wheel

| Wheel Style | Handrim Attaching Points | Single Wheel Assembly | Tire | Tube / Insert | Bearing | Bearing Spacer | Flat Free Rim Plug |
|---|--------------------------|--------------------------------|---------|---------------|-----------|----------------|--------------------------|
| Composite Pneumatic 1763BK (20x1-3/8") | 7 | 1034126 See Note(s) B, A | 1017519 | 1017520 | 40117X024 | 1015756 | N/A |
| Composite Pneumatic with Flat Free 4763BK (20x1-3/8") | 7 | 1034127 See Note(s) B, A | 1017519 | 1023322 | 40117X024 | 1015756 | 1030936 See Note(s) C |
| Composite Urethane U1763BK (20 x 1-1/4") | 7 | 1038555 See Note(s) B, A | 1037774 | N/A | 40117X024 | 1015756 | N/A |
| Composite Pneumatic 1764BK (22 x 1-3/8") | 7 | 1027149 See Note(s) A, D, B | 1018544 | 1018545 | 40117X024 | 1015756 | N/A |
| Composite Pneumatic with Flat Free 4764BK (22 x 1-3/8") | 7 | 1027150 See Note(s) A, B, D | 1018544 | 1023323 | 40117X024 | 1015756 | 1030936 See Note(s) C |
| Composite Urethane U1764BK (22 x 1-1/4") | 7 | 1032115 See Note(s) D, A | 1028250 | N/A | 40117X024 | 1015756 | N/A |

NOTE: A - Does not include handrim. Includes tire, tube/insert and spacer.
 B - When ordering 20 and 22" wheels, a wheel lock shoe spacer may also be required. See section titled, "Wheel Locks."
 C - Wheels with flat free inserts have a plug which caps an opening in the rim.
 D - 22" rear wheels are not available on 10" or 12" deep non-reclining chairs.