


Shroud Covers

R2 jr. (14" Wide)
(Before 3/1/00)

NOTE: After 10/05/99 item 5 and 6 shown were replaced with dual lock fastener strips of 2 inches long instead of 1 inch long. Right/left side panel shrouds and trunk lid were effected.



Shroud Covers

R2 jr. (14" Wide)

(Before 3/1/00)

Item #'s 1-3 Shroud Covers

R2 jr. (14" Wide)

Item Number	Description	Bright Plum Metallic	Bright Blue Metallic	Blue Green Chameleon	Candy Apple Red	Bright Silver Metallic
1	Side Panel with Rubber Bushing (Sling Seat) - Right	1087147-NLA See Note(s) A, B	1087148-NLA See Note(s) A, B	1087149-NLA See Note(s) B, A	1087150-NLA See Note(s) B, A	1087151-NLA See Note(s) A, B
1	Side Panel with Rubber Bushing (Sling Seat) - Left	1087152 See Note(s) A, B	1087153-NLA See Note(s) A	1087154-NLA See Note(s) A	1087155-NLA See Note(s) A	1087156-NLA See Note(s) A
2	Battery Cover	1085572-NLA See Note(s) B	1085573-NLA See Note(s) B	1085574-NLA See Note(s) B	1085575-NLA See Note(s) B	1086623-NLA See Note(s) B
3	Trunk Lid	1086181-NLA See Note(s) B	1086182-NLA See Note(s) B	1086183-NLA See Note(s) B	1086184-NLA See Note(s) B	1086185-NLA See Note(s) B

NOTE: A - Decals not included. See page titled, "Packaging, Literature and Decals."
B - These Shrouds are no longer available. See page titled, "Shroud Covers - After 2/29/00" for color conversion options.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1-3		See Chart	Item #'s 1-3 Shroud Covers	1	1
4		1070555-NLA	NLA - Fastener, Dual Lock (5/8" W x 2" L)	1	2
5		1070587	Washer, Flat (3/16 x 5/8 x 1/16")	1	8
6		1028755	Screw, Phillips Pan Head (#10-32 x 5/8")	1	8
7		7930332-NLA	NLA - Locknut, High Profile (#10-32)	25	1
8		1070594-NLA	NLA - Bushing, Rubber	1	2
9		1070575-NLA	NLA - Screw, SEMS w/ Washer (#10-32 x 7/8")	1	2
10		1050161-NLA	NLA - Screw, Hex Head w/ Patch (5/16-18 x 1")	1	3
11		1088261-NLA	NLA - Bracket, Shroud, Black - Right	1	1
11		1088262	Bracket, Shroud, Black - Left	1	1
12		1070577-NLA	NLA - Plate, Trunk Clip	1	2
13		1070588-NLA	NLA - Screw, Phillips Truss Head (#10-32 x 1/2")	1	4
14		7930332-NLA	NLA - Locknut, High Profile (#10-32)	25	1
15		1088497-NLA	NLA - Fastener, Panel (5/8" W x 3/4" L)	1	1
16		1070556-NLA	NLA - Fastener, Dual Lock (1" W x 2" L)	1	2
	A,B	1086186	Hardware, Side Panel Attaching	1	1

NOTE: A - Not Shown.
B - Includes Items 4-7 and 9-11.