


# Seat Frame and Hardware

(Junior 12-16" Wide)  
(Adult 14-15" Wide)



# Seat Frame and Hardware

**(Junior 12-16" Wide)  
(Adult 14-15" Wide)**

| Item Number | Note | Part Number | Description                        | Quantity Per |          |
|-------------|------|-------------|------------------------------------|--------------|----------|
|             |      |             |                                    | Package      | Assembly |
| 1           |      | 1037610     | Screw, Hex Head (1/4-20 x 1-5/8")  | 1            | 4        |
| 2           |      | 1026099     | Spacer, Coved (1/4 x 11/16 x 1/2") | 1            | 4        |
| 3           |      | 1047449-NLA | NLA - Spacer (1/4 x 3/8 x 3/16")   | 1            | 4        |
| 4           |      | 1025736     | Locknut (1/4-20)                   | 1            | 4        |
| 5           |      | 47471X074   | Spring, Arm Lock                   | 1            | 1        |
| 5           |      | 1028594     | Pin, Arm Lock, Black               | 1            | 1        |
| 5           |      | 1090768     | Kit, Lock Handle w/Rivets (Pair)   | 2            | 1        |
| 6-8         |      | See Chart   | #'s 6-8 Seat Frame and Hardware    | 1            | 1        |

## #'s 6-8 Seat Frames

| Item Number | Description                             | Junior<br>12-16" Wide<br>12-16" Deep | Adult<br>14-15" Wide<br>16-17" Deep | Adult<br>14-15" Wide<br>18-19" Deep | Adult<br>14-15" Wide<br>20-21" Deep | Adult<br>14-15" Wide<br>22" Deep |
|-------------|---|--------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|----------------------------------|
| 6           | Seat Frame<br>Right<br>See Note(s) A    | 1070473-NLA                          | 1070473-NLA                         | 1070526-NLA                         | 1070528-NLA                         | 1070530-NLA                      |
| 6           | Seat Frame<br>Left<br>See Note(s) A     | 1070474-NLA                          | 1070474-NLA                         | 1070527-NLA                         | 1070529-NLA                         | 1070531-NLA                      |
| 7           | Seat Frame<br>Center                    | 1070476-NLA                          | 1070476-NLA                         | 1070476-NLA                         | 1070533-NLA                         | 1070533-NLA                      |
| 8           | Seat Frame<br>Assembly<br>See Note(s) B | 1070469-NLA                          | 1070469-NLA                         | 1070534-NLA                         | 1070535-NLA                         | 1070536-NLA                      |

NOTE: A - Seat frame includes item 5 arm lock assembly.  
B - Includes Items 1-7