


Battery Box Tray for Group 22NF Batteries



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1092190-NLA	NLA - Battery Tray with Strap Assembly	1	1
2		1028769	Locknut (5/16-18)	1	4
3		1042484-NLA	NLA - Screw, Hex Head (5/16-18 x 2")	1	2
4		1049846	Strap Assembly, Battery	1	1
5		1074133	Washer, Steel (5/16 x 3/4 x 1/16")	1	2
NOTE: A - Battery tray design changed 9/13/99 and serves as a direct replacement for all chairs manufactured before 9/13/99 as well.					