


MKIV Mid-Line Swing Bracket Assembly



MKIV Mid-Line Swing Bracket Assembly

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1035528-NLA	NLA - Tube, Adjustment, Inside, Mid-Line, Black	1	1
2		1016778-NLA	NLA - Spring, Extension (9/16 x 1-3/8" L)	1	1
3		1035860-NLA	NLA - Tube, Adjustment, Outside, Mid-Line, Black	1	1
4		1027274	Adjustment Lock Assembly, Black	1	1
5	A	1035852-NLA	NLA - Swing Bracket Assembly, Mid-Line, Black	1	1
6		7930349-NLA	NLA - Screw, Hex Head (1/4-20 x 1-1/8")	25	1
7		1027274	Adjustment Lock Assembly, Black	1	1
8	B	1071879-NLA	NLA - Bracket, Joystick Mounting, Black	1	1
9		1027266	Half Clamp, Non-Threaded Holes	1	1
10		1027267-NLA	NLA - Half Clamp, Threaded Holes	1	1
11	B,C	1071880-NLA	NLA - Hardware, Attaching - Mid-Line	1	1
12		1049582-NLA	NLA - Bracket, Joystick Mounting, Black	1	1
13	D	1049591-NLA	NLA - Hardware, Joystick Tube Mounting	1	1

NOTE: A - Includes Items 1-4.
 B - Left side joystick mounting bracket or hardware must be ordered as a special.
 C - Includes Items 6-10.
 D - Includes Items 6, 7, 9, 10, and 12.