


Pivot Slide Tube with Footplate Assembly - J934A-14



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,B	1020234-NLA	NLA - Pivot/Slide Tube w/ Aluminum Footplate Assembly, Chrome - Right	1	1
1	A,B	1020235-NLA	NLA - Pivot/Slide Tube w/ Aluminum Footplate Assembly, Chrome - Left	1	1
2		44200X059	Footplate, Aluminum, Small, Black (6-3/4" W x 6" L)	1	1
2		44200X027	Footplate, Aluminum, Medium, Black (7-3/4" W x 6" L)	1	1
3		1016346	Cap, End, Tube, Black (7/8")	1	1
4		40120X001	Screw, Phillips Pan Head (#8-32 x 3/8")	1	2
5		1006981	Spring, Tension	1	1
6		1027097	Plug Button, Black (3/4")	1	2
7		7500005-NLA	NLA - Pivot/Slide Tube w/ Plug Button, Chrome	1	1
7		04470X013-NLA	NLA - Pivot/Slide Tube, Chrome	1	1

NOTE: A - 11-15" knee-to-heel range.
B - 14" chair width.