



Shrouds and Attaching Hardware Before 7/14/03



Side Shrouds



Battery Shrouds

Shrouds and Attaching Hardware

Before 7/14/03

#1 Side Shrouds

(Includes Items 7 and 8)

(Before 8/29/02) - Metallic Black

Right	Left
1109205-NLA	1109206-NLA
NOTE: These shrouds are no longer available. Please refer to the Magic Black chart for pair no. 1112825 as replacements.	

(After 8/28/02) - Magic Black

Right	Left	Pair
1112817-NLA	1112818-NLA	1112825-NLA
See Note(s) A	See Note(s) A	
NOTE: A - No longer available. Use 1112825 as its replacement.		

(Before 7/24/02) - Dark Blue Metallic

Right	Left
1109207-NLA	1109208-NLA
NOTE: These shrouds are no longer available. Please refer to the Dark Blue Pearl chart for pair no. 1112823 as replacements.	

(After 7/23/02) - Dark Blue Pearl

Right	Left	Pair
1112811-NLA	See Note	1112823
See Note(s) A		
NOTE: A - No longer available as a service part in the field.		

(Before 8/9/02) - Candy Apple Red

Right	Left
1109203-NLA	1109204-NLA
NOTE: These shrouds are no longer available. Please refer to the Candy Red Pearl chart for pair no. 1112824 as replacements.	

Shrouds and Attaching Hardware Before 7/14/03

(After 8/8/02) - Candy Red Pearl

Right	Left	Pair
1112805-NLA	1112806-NLA	1112824-NLA

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	#1 Side Shrouds Candy Apple Red (Before 8/9/02)	1	1
1		See Chart	#1 Side Shrouds Candy Red Pearl (After 8/8/02)	1	1
1		See Chart	#1 Side Shrouds Dark Blue Metallic (Before 7/24/02)	1	1
1		See Chart	#1 Side Shrouds Dark Blue Pearl (After 7/23/02)	1	1
1		See Chart	#1 Side Shrouds Magic Black (After 8/28/02)	1	1
1		See Chart	#1 Side Shrouds Metallic Black (Before 8/29/02)	1	1
2		1025451	Screw, Slotted Hex Washer Head (#10-32 x 5/8")	1	8
3		1070587	Washer, Flat (3/16 x 5/8 x 1/16")	1	12
4		1028755	Screw, Phillips Pan Head (#10-32 x 5/8")	1	4
5		1106366-NLA	NLA - Plate, Shroud -Right	1	1
5		1106367-NLA	NLA - Plate, Shroud -Left	1	1
6	B	1106343-NLA	NLA - Hardware, Shroud, Side (1 Pair)	1	1
7		1109032-NLA	NLA - Decal, Shroud - Pronto M71	1	1
8		1104797	Decal, Drive/Push	1	1
9	C,E	1108464-NLA	NLA - Kit, Shroud Clamp Hardware - Right	1	1
9	C,E	1108465-NLA	NLA - Kit, Shroud Clamp Hardware - Left	1	1
10		1109036-NLA	NLA - Bracket, Clamp Shroud Mounting, Black - Right	1	1
10		1109037-NLA	NLA - Bracket, Clamp Shroud Mounting, Black - Left	1	1
11		1026907	Washer (1/4 x 1/2 x 1/16")	1	1
12		1025736	Locknut (1/4-20)	1	1
13		1109035-NLA	NLA - Clamp, Shroud	1	1
14		1025713-NLA	NLA - Screw, Hex Head (1/4-20 x 1")	1	1
15	D	1106503-NLA	NLA - Battery Shroud without Decal, Black - Rear	1	1
16		1119281-NLA	NLA - Plate, Name, Patent Numbers	1	1
17	A	1106502-NLA	NLA - Battery Shroud without Decal, Black - Front	1	1
18		1095489	Decal, Invacare 3-D Logo (2-1/2" L)	1	2
19		1070587	Washer, Flat (3/16 x 5/8 x 1/16")	1	2
20		1117240	Screw, Thumb (#10-32 x 1/2")	1	2

NOTE: A - This shroud encompasses the controller as well as the batteries.
 B - Includes items 2-5
 C - Includes items 10-14
 D - This shroud encompasses the battery charger as well as the batteries.
 E - This Washer is not required for use after 2/6/07; Kit will not contain the Washer after 2/6/07.