


Solid Folding Seats



Solid Folding Seats

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1016771-NLA	NLA - Hook, Seat	1	2
2		1016770-NLA	NLA - Bracket, Seat, Chrome	1	2
3		1013218	Washer (#10 Flat)	1	8
4		7930360-NLA	NLA - Package, Screw, Phillips Pan Head (#10-32 x 11/16")	25	8
5		See Chart	# 5 Solid Folding Seats	1	1

Solid Folding Seat with Hardware

Part Number	Item #5 Solid Folding Seat	Quantity
8881036531-NLA See Note(s) A	10" Wide x 10" Deep	1
8881036532-NLA See Note(s) A	10" Wide x 11" Deep	1
8881036533-NLA See Note(s) A	10" Wide x 12" Deep	1
8881036536-NLA See Note(s) A	12" Wide x 12" Deep	1
8881036537-NLA See Note(s) A	12" Wide x 13" Deep	1
8881036538-NLA See Note(s) A	12" Wide x 14" Deep	1
8881036541-NLA See Note(s) A	14" Wide x 12" Deep	1
8881036542-NLA See Note(s) A	14" Wide x 13" Deep	1
8881036543-NLA See Note(s) A	14" Wide x 14" Deep	1

NOTE: A - For upholstery color, see Upholstery Color Page.