


Amputee Attachment



Amputee Attachment

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	47183X037-NLA	NLA - Spacer (1/2 x 7/8 x 7/8")	1	1
2		1020512-NLA	NLA - Bracket, Amputee - Right	1	1
2		1020513-NLA	NLA - Bracket, Amputee - Left	1	1
3		7930352-NLA	NLA Locknut (7/16-20)	25	1
4		7930364-NLA	NLA - Package, Locknut (1/4-20)	25	2
5		1007269	Screw, Socket Button Head (1/4-20 x 5/8")	1	2
6		1016289	Washer, Caplugs (15/32" I.D. x 1-3/32" O.D.)	1	1
7		47203X025	Bolt, Hex Head Cap (7/16-20 x 4-3/4")	1	1
8		1016291	Cap, Hub, Black	1	1
9		1016980	Spacer, Wheel Lock Shoe (2-63/64")	1	2
10		1008438	Screw, Set Socket Cup Point (#8-32 x 3/16")	1	2
11	1036618-NLA	NLA - Attachment, Amputee (Pair)	2	1	
NOTE: A - For 20" Rear Wheel					