


Torso Support



Torso Support

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		60284X001	Screw, Socket Head (1/4-20 x 1/2")	1	2
2		03018X043	Bracket, Clamp, Black- Right	1	1
2		03018X046	Bracket, Clamp, Black - Left	1	1
3		03018X012-NLA	NLA - Rod, Support w/ Guide Tube	1	1
4		7930360-NLA	NLA - Package, Screw, Phillips Pan Head (#10-32 x 11/16")	25	2
5	B,C	88804529X040-NLA	NLA - Pad, Headrest	1	1
6		03018X016-NLA	NLA - Bracket, Support, Chrome	1	1
7	B	88804529X043-NLA	NLA - Cover Torso Bracket	--	--
8	B,D	8881050189-NL A	NLA - Torso Support Assembly (Pair)	2	1
	A,E	1023610-NLA	NLA - Instruction Sheet	1	1
NOTE: A - Not Shown B - Specify upholstery color, see page titled "Frame Finishes and Upholstery Colors" C - Vinyl Only D - Non-Recliner E - Item is no longer available from the Invacare Service Parts Department. Refer to the Invacare Website under "Product Literature Search". Enter the part number in the "Keywords" field and download a free PDF version copy.					