


# Base Frame and Van Seat Post Assembly

## Before 7/14/03



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		40117X006	Bearing (5/8")	1	4
2		1105026	Plug Button, Head Tube, Black	1	2
3		1109484-NLA	NLA - Post, Adjustable Seat	1	1
4		60355X000-NLA	NLA - Locknut (3/8-24)	1	1
5	A	1110603-NLA	NLA - Base Frame Assembly	1	1
6		1111704-NLA	NLA - Screw, Socket Head (3/8-24 x 2-1/4")	1	1
7		1110539-NLA	NLA - Decal, Warning, Battery Cable	1	1
8	B	1112268-NLA	NLA - Base Frame and Adjustable Seat Post Assembly	1	1

NOTE: A - Includes items 1 and 7  
B - Includes items 2-7