



# Shrouds and Attaching Hardware

## Before 7/14/03



**Side Shrouds**



**Front and Rear Battery Shrouds**

# Shrouds and Attaching Hardware

Before 7/14/03

**Item #1 Side Shrouds  
(Includes Items 7 and 8)**

(Before 7/24/02)

**Dark Blue Metallic**

ProntoM50 Right	ProntoM50 Left	ProntoM51 Right	ProntoM51 Left
1110618-NLA	1110619-NLA	1110621-NLA	1110622-NLA

Note: When side shrouds become no longer available, refer to the side shroud pair number listed after 7/23/02.

(After 7/23/02)

**Dark Blue Pearl**

ProntoM50 Right	ProntoM50 Left	ProntoM50 Pair	ProntoM51 Right	ProntoM51 Left	ProntoM51 Pair
1112813-NLA	1112814-NLA	1112819-NLA	1112815-NLA	1112816-NLA	1112821-NLA See Note(s) A

NOTE: A - No longer available. Use either 1112815 or 1112816 as the replacement.

(Before 8/9/02)

**Candy Apple Red**

ProntoM50 Right	ProntoM50 Left	ProntoM51 Right	ProntoM51 Left
1110606-NLA	1110607-NLA	1110608-NLA	1110609-NLA

Note: When side shrouds become no longer available, refer to the side shroud pair number listed after 8/8/02.

(After 8/8/02)

**Candy Red Pearl**

ProntoM50 Right	ProntoM50 Left	ProntoM50 Pair	ProntoM51 Right	ProntoM51 Left	ProntoM51 Pair
1112807-NLA	1112808-NLA	1112820-NLA	1112809-NLA	1112806-NLA	1112822-NLA

## Shrouds and Attaching Hardware Before 7/14/03

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	#1 Side Shrouds Dark Metallic Blue Before 7/24/02	1	1
1		See Chart	#1 Side Shrouds Candy Red Pearl After 8/8/02	1	1
1		See Chart	#1 Side Shrouds Dark Blue Pearl After 7/23/02	1	1
1		See Chart	#1 Side Shrouds Dark Metallic Blue Before 7/24/02	1	1
2		1025451	Screw, Slotted Hex Washer Head (#10-32 x 5/8")	1	8
3		1070587	Washer, Flat (3/16 x 5/8 x 1/16")	1	12
4		1028755	Screw, Phillips Pan Head (#10-32 x 5/8")	1	4
5		1106366-NLA	NLA - Plate, Shroud -Right	1	1
5		1106367-NLA	NLA - Plate, Shroud -Left	1	1
6	A	1106343-NLA	NLA - Hardware, Shroud, Side (1 Pair)	2	1
7		1110604-NLA	NLA - Decal, Shroud - M50	1	1
7		1110605	Decal, Shroud - M51	1	1
8		1104797	Decal, Drive/Push	1	1
9	B	1108464-NLA	NLA - Kit, Shroud Clamp Hardware - Right	1	1
9	B	1108465-NLA	NLA - Kit, Shroud Clamp Hardware - Left	1	1
10		1109036-NLA	NLA - Bracket, Clamp Shroud Mounting, Black - Right	1	1
10		1109037-NLA	NLA - Bracket, Clamp Shroud Mounting, Black - Left	1	1
11		1026907	Washer (1/4 x 1/2 x 1/16")	1	1
12		1025736	Locknut (1/4-20)	1	1
13		1109035-NLA	NLA - Clamp, Shroud	1	1
14		1025713-NLA	NLA - Screw, Hex Head (1/4-20 x 1")	1	1
15		1095489	Decal, Invacare 3-D Logo (2-1/2" L)	1	2
16		1106502-NLA	NLA - Battery Shroud without Decal, Black - Front	1	1
16		1109480-NLA	NLA - Battery Shroud without Decal, Black - Rear	1	1
17		1109033-NLA	NLA - Knob, Shroud, Black	1	1
18		1109034-NLA	NLA - Plate, Latch Shroud	1	1
19		1108475-NLA	NLA - Washer, Wave (17/32 x 3/4 x 1/64")	1	1
20		1109031-NLA	NLA - Washer, Steel (13/64 x 1 x 3/64")	1	1
21		1113438-NLA	NLA - 3Screw, Socket Head (#10-32 x 1/2")	1	1
22	C	1112284-NLA	NLA - Hardware, Shroud Mounting - Rear	1	1
23		1070587	Washer, Flat (3/16 x 5/8 x 1/16")	1	2
24		1028755	Screw, Phillips Pan Head (#10-32 x 5/8")	1	2

NOTE: A - Includes items 2-5  
B - Includes items 10-14  
C - Includes items 17-21