

## Miscellaneous Accessories

**THE PRODUCTION OF THE P725 IS CURRENTLY SUSPENDED. CONTACT CUSTOMER SERVICE FOR ORDERING OPTIONS.**

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
	A,F,I	P725	Kit, Crutch/Cane Holder	1	1
			Holder, Crutch/Cane - P725	1	1
			Cup, Crutch/Cane Holder, Black	1	1
			Strap, Hook and Loop - P725	1	1
			Plug Button, Black (7/8")	1	1
			Screw w/ Washer (#10 x 7/8")	1	2
	F	1035917-NLA	NLA - Heel Loop (Pair)	2	1
	F	1017621	Heel Loop - Right/Left	1	1
	B,C,F	1035907	Kit, Heel Loop Attaching Hardware (Pair)	2	1
	D,E,F	1035931	Kit, Heel Loop Attaching Hardware (Pair)	2	1
	F	1114911-NLA	NLA - Knob, Adjustment, Black (5/16-18 x 1-1/2")	1	1
	F,H	1034622-NLA	NLA - Strap, Positioning, Seat/Chest w/ Auto Style Push Button Buckle (Eyelet to Eyelet Length 69")	1	1
	F,G	1034625	Strap, Positioning, Seat/Chest w/ Auto Style Push Button Buckle (Eyelet to Eyelet Length 61")	1	1
	F	40176X010-NLA	NLA - Pump, Tire	1	1

NOTE: A - Tighten knob assembly part number 1114911 comes standard on-chair for all models  
 B - 1650 Aluminum Footplate  
 C - Consists of: 2 ea. Screw, Phillips Head (1/4-20 x 2-1/4") 2 ea. Spacer (1/4 x 11/32 x 1-3/4") 2 ea. Locknut (1/4-20)  
 D - Consists of: 2 ea. Screw, Phillips Head (1/4-20 x 2-1/4") 2 ea. Spacer (1/4 x 11/32 x 1-3/4") 2 ea. Hex Nut w/ Nylon Cap, and 2 ea. Flat Washer (1/4 x 3/4 x 3/32").  
 E - 1651 Composite Footplate  
 F - Not Shown  
 G - 18" Wide Van Seat  
 H - 20" Wide Van Seat  
 I - THE PRODUCTION OF THE P725 IS CURRENTLY SUSPENDED. CONTACT CUSTOMER SERVICE FOR ORDERING OPTIONS.