


I.V. Rod and Bracket - 1318 DLX



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1042416	Rod, I.V., Fixed Height - Right Mount	1	1
1	A	1042417-NLA	NLA - Rod, I.V., Fixed Height - Left Mount	1	1
2		04701X041	Rod, I.V., Recliner, Black	1	1
3		04701X016-NLA	NLA - Mount, Socket, Black- Right	1	1
3		04701X017-NLA	NLA - Mount, Socket, Black - Left	1	1
	B,C	21243X000-NLA	NLA - Screw, Phillips Pan Head (#6-20 x 3/8")	1	1
NOTE: A - Includes items 2 and 3. B - Not Shown. C - To be used as mounting hardware.					