



# Crutch and Cane Carrier (A1569)



## Crutch and Cane Carrier (A1569)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1032224	Carrier, Crutch and Cane	1	1
2	B	1032223	Carrier, Crutch and Cane	1	1
3		1032222	Strap, Black	1	1
4		1032219	Holder, Crutch	1	1
5		1065362	Locknut (1/4-28)	1	1
6		1026907	Washer (1/4 x 1/2 x 1/16")	1	2
7		1026099	Spacer, Coved (1/4 x 11/16 x 1/2")	1	2
8		1055349	Screw, Hex Head (1/4-28 x 2")	1	1

NOTE: A - For Square Rear Frame. Includes Items 3 and 4.

B - For Curved Rear Frame. Includes Items 3 - 8.