


I.V. Rod - Telescoping - T1322

Not for use on Super Hemi Chairs



I.V. Rod - Telescoping - T1322

Not for use on Super Hemi Chairs

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	B	1102945-NLA	NLA - Rod, I.V., Telescoping	1	1
2		1042127	I.V. Rod and Bracket Assembly	1	1
3		1102765	Plate, Adapter	1	2
4		11239X000-NLA	NLA - Screw, Hex Head (5/16-24 x 1-3/4")	1	2
5		47324X011	Bushing, Pivot (5/16 x 7/8 x 3/4")	1	1
6		1102931	Spacer, Bracket (5/16 x 7/8 x 7/8")	1	1
7		60264X001	Locknut (5/16-24)	1	1
	A	1104899-NLA	NLA - Instruction Sheet, I.V. Rod	1	1
NOTE: A - Not Shown. B - Includes Items 2-7.					