



Front Casters and Mounting Hardware



*fork shown
for reference only*



Front Casters and Mounting Hardware

Wheel Style	Composite Urethane (6 x 1)	Composite Solid Rubber (8 x 1")	Composite Pneumatic (8 x 1-1/4")	Composite Pneumatic with Flat Free Insert (8 x 1-1/4")
Single Wheel Assembly	1026956 See Note(s) A	1058137 See Note(s) A	1009481 See Note(s) A	1030773 See Note(s) A
Tire	--	1040144	40176X040	40176X040
Tube/Insert	--	--	40176X041	1010904-NLA
Bearing	40117X026	1006049	1006049	1006049
Hub Spacer	1026172	1018982	1020592	1020592
Front Axle	7930362	57667X007-NLA	57667X007-NLA	57667X007-NLA
Axle Locknut	7930363-NLA	7930352-NLA	7930352-NLA	7930352-NLA
Axle Washer	NLA See Note(s) B	1026101	7930312	7930312
Axle Spacer	1053202 See Note(s) B	--	--	--
NOTE: A - Includes tire, tube/insert, bearing and spacer. B - As of 1/95, axle spacer (1053202) replaced both washer (16308X000) and spacer (47183X047) on 6" casters.				