


Seat Height Adjustment Slide Tubes



Seat Height Adjustment Slide Tubes

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,B	1110707-NLA	NLA - Kit, Seat Height Adjustment Slide Tubes and Hardware	1	1
2			NLA - Slide Tube Assembly, Seat Height Adjustment, Rear	1	2
3			Screw, Socket Head (5/16-18 x 1-1/2")	1	4
4			Locknut (5/16-18)	1	4
5			NLA - Slide Tube Assembly, Seat Height Adjustment, Front	1	2
1	A,C	1116615-NLA	NLA - Kit, Front/Rear Slide Tube w/ Hardware, Seat Height Adjustment (Pair)	1	1
2	D	1109019-NLA	NLA - Slide Tube, Seat Height Adjustment - Rear	1	2
2		1152529-NLA	NLA - Slide Tube Assembly, Seat Height Adjustment, Rear	1	2
3		1021344	Screw, Socket Head (5/16-18 x 1-1/2")	1	4
4		1026158	Locknut (5/16-18)	1	4
5		1109215-NLA	NLA - Slide Tube Assembly, Seat Height Adjustment, Front	1	2
5	C	1116614-NLA	NLA - Slide Tube Assembly, Seat Height Adjustment, Front	1	2

NOTE: A - Includes Items 2-5
 B - 300 Lb. capacity
 C - 400 Lb. capacity
 D - When replacing a "Slide Tube, Seat Height Adjustment - Rear" part # 1109019 with the part # 1152529 both sides must be replaced simultaneously.