

## Power Articulating Legrest Kits POWELR18 (Before 9/15/01)

**NOTE: 2G Tilt Only Systems originally ordered without a T.R.C.M.2G controller must include a T.R.C.M.2G controller and 4-way switch with hardware. Refer to the pages titled "Electronics" for MKIV T.R.C.M. 2G and MKIV T.A.C. of this book for a T.R.C.M. upgrade kit. Therefore, add the power leg kit dated after 7/26/04. See page titled "Power Articulating Legrests Kits After 7/26/04". Power Articulating Legs use sector style rigging hangers only and do not come packaged with the kits.**

### For 2GT, 2GR, 2GT/R Seating Systems When a T.R.C.M. 2G Is Not Present

#### POWELR18 Before 11/01/00

Quantity	Part Number	Description
1	1101399-NLA Kit Part Number See Note(s) B	Kit, POWELR 18
1	1099704-NLA See Note(s) B	NLA - Legrest, Power Articulating, Sector, Black, Short (13"-18" Range) - Right
1	1099705-NLA See Note(s) B	NLA - Legrest, Power Articulating, Sector, Black, Short (13"-18" Range) - Left
1	1101520-NLA	NLA - Kit, Power Leg Harness and Mounting Hardware
1	1098337-NLA	NLA - Instruction Sheet, Power Leg Installation
1	1104243-NLA See Note(s) A	NLA - Controller, MKIV T.R.C.M. 2G
1	1098341-NLA See Note(s) C	Instruction Sheet, MKIV TRCM 2G Upgrade

**NOTE:** A - After 10/31/00 software changed from 2.1 to 2.2 on the T.R.C.M. 2G controller to accommodate power articulating legrest. Seating Systems ordered before 11/01/00 requires the updated T.R.C.M. 2G controller.  
 B - Power articulating legrest and KITS before 9/15/01 are NLA. See power articulating legrest or KITS dated after 7/26/04. Legrest must be ordered in pairs to match.  
 C - Item is no longer available from the Invacare Service Parts Department. Refer to the Invacare Website under "Product Literature Search". Enter the part number in the "Keywords" field and download a free PDF version copy.

### For 2GT, 2GR, 2GT/R Seating Systems When a T.R.C.M. 2G Is Not Present

#### POWELR20 Before 11/01/00

Quantity	Part Number	Description
1	1101521-NLA Kit Part Number See Note(s) B	Kit, POWELR20
1	1099706-NLA See Note(s) B	Legrest Assembly, Power Articulating, Sector, Black, Long (15"-20" Range) - Right
1	1099707-NLA See Note(s) B	Legrest Assembly, Power Articulating, Sector, Black, Long (15"-20" Range) - Left
1	1101520-NLA	NLA - Kit, Power Leg Harness and Mounting Hardware
1	1098337-NLA	NLA - Instruction Sheet, Power Leg Installation
1	1104243-NLA See Note(s) A	Controller, MKIV T.R.C.M. 2G

# Power Articulating Legrest Kits

## POWELR18 (Before 9/15/01)

**NOTE: 2G Tilt Only Systems originally ordered without a T.R.C.M.2G controller must include a T.R.C.M.2G controller and 4-way switch with hardware. Refer to the pages titled "Electronics" for MKIV T.R.C.M. 2G and MKIV T.A.C. of this book for a T.R.C.M. upgrade kit. Therefore, add the power leg kit dated after 7/26/04. See page titled "Power Articulating Legrests Kits After 7/26/04". Power Articulating Legs use sector style rigging hangers only and do not come packaged with the kits.**

**For 2GT, 2GR, 2GT/R Seating Systems When a T.R.C.M. 2G Is Not Present**

### POWELR20 Before 11/01/00

Quantity	Part Number	Description
1	1098341-NLA See Note(s) B	Instruction Sheet, MKIV TRCM 2G Upgrade
<p><b>NOTE:</b> A - After 10/31/00 software changed from 2.1 to 2.2 on the T.R.C.M. 2G controller to accommodate power articulating legrest. Seating Systems ordered before before 11/01/00 requires the updated T.R.C.M. 2G controller.            B - Power articulating legrest and KITS before 9/15/01 are NLA. See power articulating legrest or KITS dated after 7/26/04. Legrest must be ordered in pairs to match.</p>		

# Power Articulating Legrest Kits

## POWELR18 (Before 9/15/01)

**NOTE: 2G Tilt Only Systems originally ordered without a T.R.C.M.2G controller must include a T.R.C.M.2G controller and 4-way switch with hardware. Refer to the pages titled "Electronics" for MKIV T.R.C.M. 2G and MKIV T.A.C. of this book for a T.R.C.M. upgrade kit. Therefore, add the power leg kit dated after 7/26/04. See page titled "Power Articulating Legrests Kits After 7/26/04". Power Articulating Legs use sector style rigging hangers only and do not come packaged with the kits.**

### For 2GT, 2GR, 2GT/R Seating Systems When a T.R.C.M. 2G Is Present

#### POWELR18 After 10/31/00

Quantity	Part Number	Description
1	1099952-NLA	Kit, POWELR18
1	1099704-NLA <i>See Note(s) A</i>	NLA - Legrest Assembly, Power Articulating, Sector, Black, Long (15"-20" Range) - Right
1	1099705-NLA <i>See Note(s) A</i>	NLA - Legrest Assembly, Power Articulating, Sector, Black, Long (15"-20" Range) - Left
1	1101520-NLA	NLA - Kit, Power Leg Harness and Mounting Hardware
1	1098337-NLA	NLA - Instruction Sheet, Power Leg Installation

**NOTE:** A - Power articulating legrest and KITS before 9/15/01 are NLA. See power articulating legrest or KITS dated after 7/26/04. Legrest must be ordered in pairs to match.

### For 2GT, 2GR, 2GT/R Seating Systems When a T.R.C.M. 2G Is Present

#### POWELR20 After 10/31/00

Quantity	Part Number	Description
1	1101117-NLA <i>See Note(s) A</i>	Kit, POWELR20
1	1099706-NLA <i>See Note(s) A</i>	NLA - Legrest Assembly, Power Articulating, Sector, Black, Long (15"-20" Range) - Right
1	1099707-NLA <i>See Note(s) A</i>	NLA - Legrest Assembly, Power Articulating, Sector, Black, Long (15"-20" Range) - Left
1	1101520-NLA	NLA - Kit, Power Leg Harness and Mounting Hardware
1	1098337-NLA	NLA - Instruction Sheet, Power Leg Installation

**NOTE:** A - Power articulating legrest and KITS before 9/15/01 are NLA. See power articulating legrest or KITS dated after 4/30/02. Legrest must be ordered in pairs to match.