

Power Non-Articulating Legrest Kits

POWELR18/POWELR20 (9/15/01-4/30/02)

NOTE: 2G Tilt Only Systems originally ordered without a T.R.C.M.2G controller must include a T.R.C.M.2G controller and 4-way switch with hardware. Refer to the pages titled "Electronics" for MKIV T.R.C.M. 2G and MKIV T.A.C. of this book for a T.R.C.M. upgrade kit. Therefore, add the power legrest kit dated after 7/26/04. See page titled "Power Articulating Legrest Kits After 7/26/04". Power Articulating Legs use sector style rigging hangers only and do not come packaged with the kits.

For 2GT, 2GR, 2GT/R Seating Systems When a T.R.C.M. 2G Is Not Present

POWELR 18/POWELR 20 Before 11/01/00

Quantity	Part Number	Description
1	1109315-NLA Kit Part Number See Note(s) B	Kit, POWELR 18/POWELR 20
1	1109313-NLA See Note(s) B	Legrest, Power Non-Articulating, Sector, Black (13"-20" Range) - Right
1	1109314-NLA See Note(s) B	Legrest, Power Non-Articulating, Sector, Black (13"-20" Range) - Left
1	1101520-NLA	Kit, Power Leg Harness and Mounting Hardware
1	1098337-NLA	Instruction Sheet, Power Leg Installation
1	1104243-NLA See Note(s) A	Controller, MKIV T.R.C.M. 2G
1	1098341-NLA See Note(s) C	Instruction Sheet, MKIV TRCM 2G Upgrade
<p>NOTE: A - After 10/31/00 software changed from 2.1 to 2.2 on the T.R.C.M. 2G controller to accommodate power articulating legrest. Seating Systems ordered before before 11/01/00 requires the updated T.R.C.M. 2G controller. B - Power articulating legrest and KITS before 9/15/01 are NLA. See power articulating legrest or KITS dated after 7/26/04. Legrest must be ordered in pairs to match. C - Item is no longer available from the Invacare Service Parts Department. Refer to the Invacare Website under "Product Literature Search". Enter the part number in the "Keywords" field and download a free PDF version copy.</p>		

For 2GT, 2GR, 2GT/R Seating Systems When a T.R.C.M. 2G Is Present

POWELR18/POWELR20 After 10/31/00

Quantity	Part Number	Description
1	1109316-NLA Kit Part Number See Note(s) A	Kit, POWELR18/POWELR20
1	1109313-NLA See Note(s) A	NLA - Legrest, Power Non-Articulating, Sector, Black (13"-20" Range) - Right
1	1109314-NLA See Note(s) A	NLA - Legrest, Power Non-Articulating, Sector, Black (13"-20" Range) - Left
1	1101520-NLA	NLA - Kit, Power Leg Harness and Mounting Hardware
1	1098337-NLA	NLA - Instruction Sheet, Power Leg Installation
<p>NOTE: A - Power articulating legrest and KITS before 9/15/01 are NLA. See power articulating legrest or KITS dated after 7/26/04. Legrest must be ordered in pairs to match.</p>		