

Power Articulating Legrest Kits

ELRPW (5/1/02-7/26/04)

NOTE: 2G Tilt Only Systems originally ordered without a T.R.C.M.2G controller must include a T.R.C.M.2G controller and 4-way switch with hardware. Refer to the pages titled "Electronics" for MKIV T.R.C.M. 2G and MKIV T.A.C. of this book for a T.R.C.M. upgrade kit. Therefore, add the power leg kit dated after 7/26/04. See page titled "Power Articulating Legrest kits After 7/26/04". Power Articulating Legs use sector style rigging hangers only and do not come packaged with the kits.

For 2GT, 2GR, 2GT/R Seating Systems When a T.R.C.M. 2G Is Not Present

ELRPW Before 5/1/02

Quantity	Part Number	Description
1	1126096-NLA See Note(s) B, C, C1	Kit, ELRPW Power Articulating Legrest
1	1126096E-NLA	Reconditioned - Power Articulating Legrest Kit
1	1126091 See Note(s) C, C1	Legrest Assembly, Power Articulating, Sector, Black (13"-20" Range) - Right
1	1126092 See Note(s) C1, C	Legrest Assembly, Power Articulating, Sector, Black (13"-20" Range) - Left
1	1101520-NLA	Kit, Power Leg Harness and Mounting Hardware
1	1098337-NLA	Instruction Sheet, Power Leg Installation
1	1109548-NLA See Note(s) A, A1	Controller w/ Hardware, MKIV TRCM 2G
1	1109548E-NLA	Reconditioned - Controller Assembly, MKIV TRCM 2G
1	1098341-NLA See Note(s) D	Instruction Sheet, MKIV TRCM 2G Upgrade

NOTE: A - After 4/30/02 MKIV T.R.C.M. controller software changed from 2.21 to 2.3 to accommodate power articulating legrest made after 4/30/02. A new MKIV T.R.C.M. controller must be ordered when replacing power articulating or non-articulating legrest made before
A1 - 5/1/02. Verify software version before ordering.
B - Kit includes the 6 part numbers listed below the kit part number
C - Due to a product enhancement that provides a higher pivot point at the knee, kit part number 1111723 and right/left side legrest 1110743 and 1110744 were replaced with the current part numbers shown on after 7/26/04.
C1 - Right and left side legrest must be ordered in pairs prior to 7/27/04.
D - Item is no longer available from the Invacare Service Parts Department. Refer to the Invacare Website under "Product Literature Search". Enter the part number in the "Keywords" field and download a free PDF version copy.

Power Articulating Legrest Kits

ELRPW (5/1/02-7/26/04)

NOTE: 2G Tilt Only Systems originally ordered without a T.R.C.M.2G controller must include a T.R.C.M.2G controller and 4-way switch with hardware. Refer to the pages titled "Electronics" for MKIV T.R.C.M. 2G and MKIV T.A.C. of this book for a T.R.C.M. upgrade kit. Therefore, add the power leg kit dated after 7/26/04. See page titled "Power Articulating Legrest kits After 7/26/04". Power Articulating Legs use sector style rigging hangers only and do not come packaged with the kits.

For 2GT, 2GR, 2GT/R Seating Systems When a T.R.C.M. 2G Is Present

ELRPW 5/1/02-7/26/04

Quantity	Part Number	Description
1	1111722-NLA Kit Part Number See Note(s) B, A, B1	Kit, ELRPW Power Articulating Legrest Kit
1	1110743-NLA See Note(s) B1, B	NLA - Legrest Assembly, Power Articulating, Sector, Black (13"-20" Range) - Right
1	1110744-NLA See Note(s) B, B1	NLA - Legrest Assembly, Power Articulating, Sector, Black (13"-20" Range) - Left
1	1101520-NLA	NLA - Kit, Power Leg Harness and Mounting Hardware
1	1098337-NLA	NLA - Instruction Sheet, Power Leg Installation
<p>NOTE: A - Kit includes the 4 part numbers listed below the kit part number B - Due to a product enhancement that provides a higher pivot point at the knee, a customer must order a complete right and left side legrest assembly. A kit is required when ordering Power Articulating Legrest for the first time. B1 - See page titled "Power Articulating Legrest Kits After 7/26/04".</p>		