



Pivot Slide Tube with Footplate

H93E / H93M / H93M WITH FPP - Footrests



Pivot Slide Tube with Footplate

H93E / H93M / H93M WITH FPP - Footrests

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	H	See Chart	Pivot Slide Tube with Footplate	1	1
2	A	1041528-NLA	NLA - Pivot/Slide Tube w/ Plug Button, Chrome	1	1
2	B	7500084-NLA	NLA - Pivot/Slide Tube w/ Plug Button, Chrome	1	1
2	C	1041776-NLA	NLA - Pivot/Slide Tube w/ Plug Buttons, Chrome	1	1
2	D	1031691-NLA	NLA - Pivot/Slide Tube w/ Plug Button, Chrome	1	1
3		1027097	Plug Button, Black (3/4")	1	2
4	A,C	44200X059	Footplate, Aluminum, Small, Black (6-3/4" W x 6" L)	1	1
4	B,D	1043659-NLA	NLA - Footplate, Aluminum, Deluxe, Black - Right (8" W x 5-3/4" L)	1	1
4	B,D	1043660-NLA	NLA - Footplate, Aluminum, Deluxe, Black - Left (8" W x 5-3/4" L)	1	1
5	E	1025729	Screw, Phillips Pan Head Tap (#8-32 x 3/8")	1	2
6	F	1006981	Spring, Tension	1	1
6	G	1036120-NLA	NLA - Spring	1	2
7		1027097	Plug Button, Black (3/4")	1	1
NOTE: A - For use on H93E/H93M models with a 14 - 15" chair width. B - For use on H93E/H93M models with a 16 - 22" chair width. C - For use on H93M with FPP with a 14 - 15" chair width. D - For use on H93M with FPP with a 16 - 22" chair width. E - Chairs 14-15" wide require (2) screws and chairs 16 - 22" wide require (1) screw. F - For use on chairs 14 - 15" wide. G - For use on chairs 16-22" wide. H - Includes items 2-7.					
NOTE: Pivot Slide Tube and Footplate is attached to Footrest Support using a compression clamp welded to the Support Assembly.					

Pivot Slide Tube with Footplate

Prior Model Number	Current Model Number	Chair Width	Side	Part Number
H93E, H93M	H93E, H93M - Chrome	14 - 15"	Right	1020232-NLA
H93E, H93M	H93E, H93M - Chrome	14 - 15"	Left	1020233-NLA
H93E, H93M	H93E, H93M - Chrome	16 - 22"	Right	1028905-NLA
H93E, H93M	H93E, H93M - Chrome	16 - 22"	Left	1028906-NLA
H93M with FPP	H93M with FPP - Chrome	14 - 15"	Right	1031936-NLA
H93M with FPP	H93M with FPP - Chrome	14 - 15"	Left	1031937-NLA
H93M with FPP	H93M with FPP - Chrome	16 - 22"	Right	1033381-NLA
H93M with FPP	H93M with FPP - Chrome	16 - 22"	Left	1033382-NLA