


Crossbraces and Hardware



Crossbraces and Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	#1 Crossbraces without Hardware	1	1
2	A	1100791	Hardware, Crossbrace (20" W)	1	1
2	A	1100828	Hardware, Crossbrace (22" W)	1	1
2	A	1100829	Hardware, Crossbrace (24" W)	1	1
2	A	1100830	Hardware, Crossbrace (26" W)	1	1
2	A	1100831	Hardware, Crossbrace (28" W)	1	1
2	A	1100832	Hardware, Crossbrace (30" W)	1	1
3		1027095	Plug Button, Black (7/8")	1	4
4		1008495	Screw, Socket Button Head (1/4-20 x 1-1/2")	1	10
5		1100833	Link, Black (20" W)	1	2
5		1100834	Link, Black (22" W)	1	2
5		1100835	Link, Black (24" W)	1	2
5		1100836	Link, Black (26" W)	1	2
5		1100837	Link, Black (28" W)	1	2
5		1100838-NLA	NLA - Link, Black (30" W)	1	2
6		60264X001	Locknut (5/16-24)	1	1
7		40317X001	Screw, Hex Head (5/16-24 x 2-3/4")	1	2
8		1012180	Bracket, Pivot, Black	1	2
9		60275X000	Locknut (1/4-20)	1	1
10		1012172-NLA	NLA - Clamp, Saddle	1	4

NOTE: A - Includes Items 3-10.

#1 Crossbraces without Hardware

Chair Width:	Side	18" Chair Depth	20" Chair Depth
20"	Right	1098803	1098805
	Left	1098804	1098806
22"	Right	1098834	1098836
	Left	1098835	1098837
24"	Right	1098839	1098841
	Left	1098840	1098842
26"	Right	1098844	1098846
	Left	1098845	1098847
28"	Right	1098849	1098851
	Left	1098850	1098852
30"	Right	1098782	1098790
	Left	1098783	1098791