


Back Canes and Hardware



Back Canes and Hardware

Back Canes and Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1100785	Kit, Back Cane Hardware (Pair)	1	1
2			Plug Button, Black (7/8")	1	2
3			Grip, Foam Hand (6-7/8" L)	1	2
4			Locknut (5/16-18)	1	4
5			Screw, Socket Head (5/16-18 x 1-1/2")	1	4
6		1100769	Cane, Back, Black	1	1

NOTE: A - Includes items 2-5. Hardware quantity accommodates the right and left side frames