


Frame

Please call Customer Service when ordering replacement frames and frame hardware, these are special orders and must be ordered through Customer Service.
Please refer to specials number C-050922110721 when ordering.




A4 Frame Order Requirements


Frame orders must be treated as a custom order.
When placing the order, the following information is required:

1. Tapered or "V" Front.
2. Width and Depth.
3. Front seat to floor.
4. Rear seat to floor.
5. Caster size.
6. Rear wheel size.
7. Frame color.

Tapered Front End



"V" Front End



Frame

Please call Customer Service when ordering replacement frames and frame hardware, these are special orders and must be ordered through Customer Service.

Please refer to specials number C-050922110721 when ordering.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,B	8881114110-NL A	NLA - Frame, A-4 Aluminum, Color	1	1
1	A,B	8881116378-NL A	NLA - Frame, A-4 Titanium, Color	1	1
NOTE: A - See frame ordering requirements. B - Seat frame not included. See category titled, "Seat Frame and Frame Hardware."					